

# SALTER-STYLE® 16SOFT HIGH FLOW CANNULA

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The comfort and convenience of oxygen therapy with nasal cannulas are now readily available for patients requiring higher oxygen flow rates. AirLife™ presents the 16SOFT High Flow Cannula from Salter Labs® that can replace cumbersome face masks with its enhanced comfort facepiece and 3-Channel safety tubing for flow rates up to 15 L/min.



## FEATURES

- Delivers O<sub>2</sub> flows up to 15 L/min with 3-Channel safety tubing
- Patented slide bolo
- Enhanced reservoir facepiece
- 16SOFT High Flow Cannulas available in 7', 14' and 25' tubing lengths

## PROBLEMS SOLVED

- Eliminates the need to use face masks
- Increased patient mobility
- Less cumbersome and claustrophobic than a mask
- Simple, secure adjustments

# SALTER-STYLE® HIGH FLOW CANNULA

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

High Flow 1600SOFT-HF Cannulas are anatomically designed and shaped with soft, curved nasal prongs to provide long-term comfort.

## ADULT HIGH FLOW CANNULAS (FOR OXYGEN FLOWS UP TO 15 L/MIN)

ITEM	DESCRIPTION	TUBING LENGTH	CONNECTOR	PK
16SOFT-HF-7-25	Salter-Style Soft High Flow Cannula with Soft Tubing	7'	Standard Ribbed	25
16SOFT-HF-14-25	Salter-Style Soft High Flow Cannula with Soft Tubing	14'	Standard Ribbed	25
16SOFT-HF-25-10	Salter-Style Soft High Flow Cannula with Soft Tubing	25'	Standard Ribbed	10