

## SALTER-STYLE® SOFT OXYGEN NASAL CANNULAS

DEHP FREE | DISPOSABLE | INDIVIDUAL PACKAGING | LATEX FREE | SINGLE PATIENT MULTIPLE USE

- Comfortable, soft curved prongs gently fit into nostrils without pressing on upper lip
- High oxygen flow up to 15 LPM
- Low oxygen flow up to 3 LPM and 6 LPM
- Soft, pliable, stretchable tubing curves around cheeks and ears, lessening any pressure



### SALTER-STYLE SOFT OXYGEN NASAL CANNULAS


ITEM	DESCRIPTION	SIZE	STYLE	LENGTH	EARPIECE	O <sub>2</sub> CONNECTOR	FLOW	PK
16SOFTTG-INF-7-50	Salter-Style Soft Oxygen Cannula	Infant	Standard	7'	Soft Ears	Thread Grip	up to 3 LPM	50
16SOFT-4-50	Salter-Style Soft Oxygen Cannula	Adult	Standard	4'	Soft Ears	Trumpet	up to 6 LPM	50
16SOFT-7-50	Salter-Style Soft Oxygen Cannula	Adult	Standard	7'	Soft Ears	Trumpet	up to 6 LPM	50
16SOFTTG-7-50	Salter-Style Soft Oxygen Cannula	Adult	Standard	7'	Soft Ears	Thread Grip	up to 6 LPM	50
16SOFT-14-50	Salter-Style Soft Oxygen Cannula	Adult	Standard	14'	Soft Ears	Trumpet	up to 6 LPM	50
16SOFT-HF-7-25	Salter-Style Soft Oxygen Cannula	Adult	High Flow	7'	Soft Ears	Trumpet	up to 15 LPM	25

## SALTER-STYLE® OXYGEN NASAL CANNULAS

DISPOSABLE | INDIVIDUAL PACKAGING | LATEX FREE | SINGLE PATIENT MULTIPLE USE

The Salter-Style® Oxygen Nasal Cannulas come with a curved, soft facepiece for patient comfort.

- 3-Channel safety tubing
- Comfortable soft curved facepiece
- For relief of pressure above the ears
- Optional TLC ear wraps for relieving ear pressure



### SALTER-STYLE ADULT HIGH FLOW OXYGEN NASAL CANNULAS

ITEM	DESCRIPTION	SIZE	STYLE	TUBING LENGTH	EARPIECE	FLOW	O <sub>2</sub> CONNECTOR	PK
1600-4-50	Salter-Style Oxygen Cannula	Adult	Standard	4'	Standard	up to 6 LPM	Trumpet	50
1600-7-50	Salter-Style Oxygen Cannula	Adult	Standard	7'	Standard	up to 6 LPM	Trumpet	50
1600HF-7-25	Salter-Style Oxygen Cannula	Adult	High Flow	7'	Standard	up to 15 LPM	Trumpet	25
1606B-0-50	Salter-Style Oxygen Cannula	Adult	Standard	NA	Standard	up to 6 LPM	Barbed	50
1600TLC-7-25	Salter-Style Oxygen Cannula	Adult	Standard	7'	TLC Ear Wraps	up to 6 LPM	Trumpet	25
1601-7-50	Salter-Style Oxygen Cannula	Infant	Low Flow	7'	Standard	up to 3 LPM	Trumpet	50